



Human Noradrenaline ELISA Kit

Enzyme Immunoassay for the quantification of Noradrenaline in human plasma and urine samples.

Catalog number: ARG80472

Package: 96 wells

For research use only. Not for use in diagnostic procedures.

TABLE OF CONTENTS

SECTION	Page
INTRODUCTION	3
PRINCIPLE OF THE ASSAY	4
MATERIALS PROVIDED & STORAGE INFORMATION	5
MATERIALS REQUIRED BUT NOT PROVIDED	6
TECHNICAL HINTS AND PRECAUTIONS	6
SAMPLE COLLECTION & STORAGE INFORMATION	8
REAGENT PREPARATION	8
ASSAY PROCEDURE	9
CALCULATION OF RESULTS	11
EXAMPLE OF TYPICAL STANDARD CURVE	12
QUALITY ASSURANCE	13

MANUFACTURED BY:

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: info@arigobio.com

INTRODUCTION

In humans the catecholamines adrenaline (epinephrine), noradrenaline (norepinephrine) and dopamine are neurotransmitters of the sympathetic nervous system and are involved in many physiological processes. The sympathetic nervous system sets the body to a heightened state of alert, also called as the body's fight-or-flight response.

In the human body the catecholamines and their metabolites indicate the adaption of the body to acute and chronic stress.

Next to the metanephrine/normetanephrine the catecholamines are important for the diagnosis and the follow-up of tumors of the sympathoadrenal system like the pheochromocytomas. The quantitative determination of catecholamines in urine is preferred for the diagnosis of these tumors, whereas the determination of catecholamines in plasma is medically sensible for the localization of the tumor and for function testing. Values above the cut-off can provide an indication for neuroendocrine tumors. However, in literature various diseases like hypertension, cardiovascular diseases, schizophrenia and manic depression are described with abnormal low or high levels of catecholamines. Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as under point "Procedural cautions, guidelines and warnings". Any laboratory result is only a part of the total clinical picture of the patient. Only in cases where the laboratory results are in an acceptable agreement with the overall clinical picture of the patient it can be used for therapeutic consequences.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

PRINCIPLE OF THE ASSAY

This is an Enzyme Immunoassay for the quantification Noradrenaline in plasma and urine samples.

This assay employs the competitive quantitative enzyme immunoassay technique. Noradrenaline are first extracted by using a cis-diol-specific affinity gel, acylated and derivatized enzymatically.

The antigen has been pre-coated onto a microtiter plate. Extracted and derivatized controls, standards or samples and the solid phase bound analytes compete for a fixed number of antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum complexes are removed by washing. Anti-rabbit IgG conjugated to Peroxidase is added to each well and incubated. After washing away any unbound antibody conjugate, a substrate solution (TMB) is added to the wells and color develops in inverse-proportion to the amount of Noradrenaline present in the samples. The color development is stopped by the addition of STOP solution and the intensity of the color is measured at a wavelength of $450\text{nm} \pm 2\text{nm}$. The concentration of Noradrenaline in the sample is then determined by comparing the O.D of samples to the standard curve.

MATERIALS PROVIDED & STORAGE INFORMATION

Store the unopened kit at 2-8 °C. Use the kit before expiration date.

Component	Quantity	Storage information
Noradrenaline-coated microplate	12 strips X 8 wells	4°C
Adhesive foil	4 pieces	RT
50X Wash Buffer	20 ml	4°C
Anti-rabbit IgG-peroxidase conjugate	12 ml (Ready-to-use)	4°C
TMB substrate	12 ml (Ready-to-use)	4°C (Protect from light)
STOP solution	12 ml (Ready-to-use)	4°C
Standard A-F (0, 5, 20, 75, 250, 1000 ng/ml)	4 ml each (Ready-to-use)	4°C
Noradrenaline Antiserum	6 ml (Ready-to-use)	4°C
Adjustment Buffer	4 ml (Ready-to-use)	4°C
Acylation Buffer	20 ml (Ready-to-use)	4°C
Acylation Reagent	3 ml (Ready-to-use)	4°C
Assay Buffer	6 ml (Ready-to-use)	4°C
Coenzyme (S-adenosyl-L-methionine)	4 ml (Ready to use)	4°C
Enzyme (COMT)	2 vials (Lyophilized)	4°C
Extraction Buffer	6 ml (Ready to use)	4°C
Extraction Plate (coated with boronate affinity gel)	2 X 48 wells (Ready-to-use)	4°C
Control 1 (40 ng/ml; \pm 40%)	4 ml (Ready-to-use)	4°C
Control 2 (150 ng/ml; \pm 40%)	4 ml (Ready-to-use)	4°C
Hydrochloric Acid	20 ml (Ready-to-use)	4°C

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable of measuring absorbance at 450nm (optional: read at 620-650 nm as the reference wave length)
- Pipettes and pipette tips
- Deionized or distilled water
- Orbital microplate shaker: 3 mm (0.1118 in) 600 ± 10 rpm or 19 mm (0.75 in) 170 ± 10 rpm.
- Automated microplate washer (optional)

TECHNICAL HINTS AND PRECAUTIONS

- This kit is intended for professional use only. Users should have a thorough understanding of this protocol for the successful use of this kit. Only the test instruction provided with the kit is valid and has to be used to run the assay. Reliable performance will only be attained by strict and careful adherence to the instructions provided.
- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- Store the unopened reagents at 2- 8°C until expiration date. Do not use components beyond the expiry date indicated on the kit labels. Once opened the reagents are stable for 1 month when stored at 2 – 8 °C.
- The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch with desiccant and used in the frame provided. Microtiter strips which are removed from the frame for usage should be marked accordingly to avoid any mix-up. Once the resealable pouch has been opened, care should be taken to close it tightly with desiccant again.

Human Noradrenaline ELISA Kit ARG80472

- All kit reagents and specimens should be brought to room temperature (20 – 25 °C) and mixed gently but thoroughly before use. Avoid repeated freezing and thawing of reagents and specimens.
- If crystals are observed in the 50X Wash buffer, warm to RT or 37°C until the crystals are completely dissolved.
- Ensure complete reconstitution and dilution of reagents prior to use.
- It is highly recommended that the standards, samples and controls be assayed in duplicates.
- Once the test has been started, all steps should be completed without interruption. Make sure that the required reagents, materials and devices are prepared ready at the appropriate time.
- The controls should be included in each run and fall within established confidence limits. The confidence limits are listed in the QC-Report.
- Change pipette tips between the addition of different reagent or samples.
- The Microtiter Strips In rare cases residues of the blocking and stabilizing reagent can be seen in the wells as small, white dots or lines. These residues do not influence the quality of the product.

SAMPLE COLLECTION & STORAGE INFORMATION

Plasma - Collect plasma using EDTA as an anticoagulant. Do not use haemolytic or lipemic samples. Assay immediately (up to 6 hours at 2-8 °C), or aliquot and store samples at ≤ -20 °C (up to 6 months). Avoid repeated freeze-thaw cycles. Avoid exposure to direct sunlight.

Urine – Spontaneous or 24-hour urine, collected in a bottle containing 10-15 ml of 6M HCl. Assay immediately (up to 48 hours at 2-8 °C, and up to 24 hours at RT), for longer storage, aliquot and store samples at ≤ -20 °C (up to 6 months). Avoid repeated freeze-thaw cycles. Avoid exposure to direct sunlight.

REAGENT PREPARATION

- **1X Wash buffer:** Dilute 50X Wash buffer into distilled water to yield 1X Wash buffer. E.g. dilute the 20 ml of 50X Wash Buffer into distilled water to a final volume of 1000 ml. The diluted wash buffer can be stored at 2 – 8 °C for 1 month.
- **Enzyme solution:** Reconstitute the lyophilized “Enzyme (COMT)” with 1ml of **distilled water** and mix well. Add **0.3 ml** of **Coenzyme** followed by **0.7 ml** of **Adjustment buffer**. The total volume of Enzyme solution is 2ml.
Note: The Enzyme Solution has to be prepared freshly prior to the assay (not longer than 10- 15 minutes in advance). Discard after use!

ASSAY PROCEDURE

All materials should be equilibrated to room temperature (RT, 20 – 25 °C) and mix thoroughly before use. Standards, samples and controls should be assayed in duplicates. If a microplate shaker is not available, please optimize the shaking speed of user's shaker.

Sample Preparation, Extraction and Acylation

1. Pipette **10 µl** of **standards, controls, urine samples** and **300 µl** of **plasma samples** into the appropriate wells of the Extraction Plate.
2. Add **250 µl** of **distilled water** to wells with **standards, controls and urine samples**.
3. Add **50 µl** of **Assay Buffer** to all wells.
4. Add **50 µl** of **Extraction Buffer** to all wells.
5. Cover plate and incubate for **30 mins at RT** on a microplate shaker (~600 rpm).
6. Remove foil, discard and blot dry by tapping the inverted plate on absorbent material. **Wash** each well with **1 ml** of **1X wash buffer** and **shake** for **5 min at RT** on a microplate shaker (600 rpm). Blot dry by tapping the inverted plate on absorbent material.
7. Repeat **wash as step 6**. Discard and blot dry by tapping the inverted plate on absorbent material.
8. Add **150 µl** of **Acylation Buffer** into all wells.
9. Add **25 µl** of **Acylation Reagent** into all wells.
10. Incubate for **15 mins at RT** on a microplate shaker (~600 rpm).
11. Remove foil, discard and blot dry by tapping the inverted plate on absorbent material. **Wash** each well with **1 ml** of **1X wash buffer** and **shake**

Human Noradrenaline ELISA Kit ARG80472

- for **10 min at RT** on a microplate (~600 rpm). Blot dry by tapping the inverted plate on absorbent material.
12. Add **150 µl** of **Hydrochloric Acid** into all wells.
 13. Cover plate and **incubate for 10 mins at RT** on a microplate shaker (~600 rpm).
 14. Remove foil, **do not decant the supernatant!**
 15. Use **20 µl** for Noradrenaline assay.

Noradrenaline ELISA procedure

1. Remove excess microtiter strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal it.
2. Add **25 µl** of **Enzyme solution** (refer to REAGENT PREPARATION) into all wells of Noradrenaline-coated Microtiter Strips.
3. Add **20 µl** of the extracted standards, controls and samples into the appropriate wells of Noradrenaline-coated Microtiter Strips.
4. Incubate for **30 mins at RT** on a microplate (~600 rpm).
5. Add **50 µl** of **Noradrenaline Antiserum** into all wells.
6. Cover plate with Adhesive foil and **incubate for 2 hours at RT** on a microplate (~600 rpm).
7. Remove the foil and discard. **Aspirate each well and wash**, repeating the process 2 times for a **total 3 washes**. Wash by filling each well with **1× Wash Buffer (300 µl)** using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating, decanting or blotting against clean paper towels.
8. Add **100 µl** of **Anti-rabbit IgG-peroxidase conjugate** into wells.

Human Noradrenaline ELISA Kit ARG80472

9. Incubate for **30 mins at RT** on a microplate (~600 rpm).
10. Aspirate each well and **wash as step 7**.
11. Add **100 µl of TMB substrate solution** into each well. **Incubate for 20-30 mins at RT** with shaking (600 rpm). Avoid exposure to light.
12. Add **100 µl of Stop Solution** to each well and shake lightly to ensure homogeneous mixing.
13. **Read** the OD with a microplate reader **at 450 nm** (optional: read with a reference wavelength between 620 nm and 650 nm) **within 10 minutes**.

CALCULATION OF RESULTS

1. Calculate the average absorbance values for each set of standards, controls and samples.
2. Using log-log, semi-log or linear graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
5. arigo provides GainData®, an in-house development ELISA data calculator, for ELISA data result analysis. Please refer our GainData® website for details. (<https://www.arigobio.com/elisa-analysis>)

Human Noradrenaline ELISA Kit ARG80472

6. The concentrations of undiluted samples and controls can be read directly from the standard curve.


7. Refer to the table below for molar conversion:

	Concentration of standards					
Standard	A	B	C	D	E	F
Noradrenaline (ng/ml)	0	5	20	75	250	1000
Noradrenaline (nmol/L)	0	30	118	443	1478	5910
Conversion	Noradrenaline (ng/ml) x 5.91 = Noradrenaline (nmol/L)					

- **Urine samples and Controls:** The concentration of urine samples and controls can be read directly from the standard curve.
- **Calculate the 24h excretion for each urine sample:** $\mu\text{g}/24\text{h} = \mu\text{g}/\text{L} \times \text{L}/24\text{h}$
- **Plasma samples:** The read concentrations have to be **divided by 30**.

EXAMPLE OF TYPICAL STANDARD CURVE

The following data is for demonstration only and cannot be used in place of data generations at the time of assay.



QUALITY ASSURANCE

Sensitivity

LOD: Urine: 1.7 ng/ml; plasma: 0.04 ng/ml

LOQ: Urine: 2.5 ng/ml; plasma: 0.093 ng/ml

Measuring range

Urine: 2.5- 1000 ng/ml

Plasma: 0.093 – 33.333 ng/ml

Assay Range

5-1000 ng/ml

Specificity

Substance	Cross Reactivity (%)
	Noradrenaline
Derivatized Adrenaline	0.08
Derivatized Noradrenaline	100
Derivatized Dopamine	0.03
Metanephrine	< 0.01
Normetanephrine	0.16
3-Methoxytyramine	< 0.01
3-Methoxy-4-hydrophenylcol	< 0.01
Tyramine	< 0.01
Phenylalanine	< 0.01
Caffeinic acid	< 0.01
L-Dopa	< 0.01
Homovanillic acid	< 0.01
Tyrosine, 3-Methoxy-4-hydroxymandelic acid	< 0.01

Intra-assay and Inter-assay precision

The CV value of intra-assay precision was 13.8% (urine) and 14.3% (plasma) and inter-assay precision was 13.2-20% (urine) and 9.2-10.9% (plasma).
[Noradrenaline]

Recovery

91-113% (Noradrenaline-urine), 75-107% (Noradrenaline-plasma)

Linearity

100-127% (Noradrenaline-urine), 102-125% (Noradrenaline-plasma)