

Human Residual RNA Detection Kit

Human Residual RNA Detection Kit is designed to detect residual Human RNA in biological products during production.

Catalog number: ARG83096

Package: 100 tests

TABLE OF CONTENTS

SECTION	Page
INTRODUCTION	3
PRINCIPLE OF THE ASSAY	3
MATERIALS PROVIDED & STORAGE INFORMATION	4
MATERIALS REQUIRED BUT NOT PROVIDED	4
TECHNICAL HINTS AND PRECAUTIONS	4
REAGENT PREPARATION	5
ASSAY PROCEDURE	6

MANUFACTURED BY:

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: info@arigobio.com

INTRODUCTION

Human Residual RNA Detection Kit is designed for the quantitative detection of residual Human total RNA in various biological products to improve control quality of nucleic acid.

Human Residual RNA Detection Kit adopts the principle of the RT-PCR fluorescent probe, combining reverse transcription PCR technology and fluorescent probe method, to realize one-step quantitative detection.

PRINCIPLE OF THE ASSAY

The Human Residual Total RNA Detection Kit is a comprehensive test kit that combines RT-PCR and qPCR technologies to detect and quantify residual Human RNA.

By utilizing reverse transcription PCR (RT-PCR), the kit converts Human RNA into complementary DNA (cDNA), enabling amplification and analysis. The subsequent quantitative PCR (qPCR) step allows for the simultaneous amplification and detection of DNA using fluorescent dyes. The Human Residual Total RNA Detection Kit includes specific primers and probes designed to target Human RNA sequences, ensuring accurate quantification of residual Human RNA in various samples.

MATERIALS PROVIDED & STORAGE INFORMATION

Store the unopened kit at -20 °C. Use the kit before expiration date.

Component	Quantity	Storage information
10x Human RNA Standard	25 μl (2ng/μl)	-20°C
RNA Dilution Buffer	2 X 1 ml	-20°C
One step RT-qPCR Buffer	1 ml	-20°C
One step Enzyme Mix	100 μΙ	-20°C
Human Primer and Probe Mix	370 μΙ	-20°C
IPC Primer and Probe Mix	370 μΙ	-20°C
ROX (High)	50 μΙ	-20°C
ROX (Low)	50 μΙ	-20°C

MATERIALS REQUIRED BUT NOT PROVIDED

- PCR machine
- Pipettes and pipette tips
- DNase/RNase-Free Water
- PCR tube

TECHNICAL HINTS AND PRECAUTIONS

- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- Store the kit at -20°C at all times.
- All reagents must be kept on ice during the entire experiment.
- Once the assay has been started, all steps should be completed without

Human Residual RNA Detection Kit ARG83096

interruption.

- It is highly recommended that the standards and samples be assayed in triplicates.
- Change pipette tips between the addition of different reagent or samples.

REAGENT PREPARATION

 Standards: Dilute 10x Human RNA Standard with RNA Dilution buffer to yield a stock concentration of 200 pg/μl. The RNA Dilution Buffer serves as zero standard (0 pg/ml), and the rest of the standard 10-fold serial Dilute Human DNA standard as according to the table below:

Standard	DNA Conc.	μl of DNA Dilution buffer	μl of standard
S6	200 pg/μl	45 μl	5 μl (stock)
S5	20 pg/μl	45 μl	5 µl (S6)
S4	2 pg/μl	45 μl	5 μl (S5)
S3	0.2 pg/μl	45 μl	5 μl(S4)
S2	0.02 pg/μl	45 μl	5 μl(S3)
S1	0.002 pg/μl	45 μl	5 μl(S2)
S0	0 pg/μl	50μΙ	0μΙ

ASSAY PROCEDURE

1 Reverse Transcription PCR mix buffer:

Total	15µl (1 wells)
ROX *	0.4μΙ
Human Primer and probe mix	3.6µl
One step Enzyme Mix	1 μΙ
One step RT-qPCR Buffer	10 μΙ

- * Choose the appropriate ROX (High or Low) to match the PCR machine. If the PCR machine does not require ROX, adjust the volume with DNase/RNase-free water to obtain a final volume of 15µl.
- 2 Mix 15μ l qPCR mix buffer with 5μ l diluent standard / sample / blank in PCR tube. The final volume should be 20 μ l.
- 3 Reverse transcription: 50°C, 15 min
- 4 Initial denaturation: 95°C, 30 sec.
- 5 PCR cycle:

95°C, 10 sec; 60°C, 40 sec, for **45 cycle**, **20 μl**.