

Product datasheet

info@arigobio.com

ARG83075 Canine Parvovirus (CPV) ELISA Kit

Package: 96 wells Store at: 4°C

Summary

Product Description ARG83075 Canine Parvovirus (CPV) ELISA Kit is an Enzyme Immunoassay kit for the qualitative of Canine

Parvovirus in dog serum.

Tested Reactivity Dog

Tested Application ELISA

Target Name Canine Parvovirus

Conjugation HRP

Conjugation Note Substrate: TMB and read at 450 nm.

Detection Range Cut - off

Sample Type Serum
Sample Volume 100 µl

Precision Intra-Assay CV: less than 7 %

Inter-Assay CV: less than 7 %

Alternate Names Canine Parvovirus, CPV, CPV2, parvo

Application Instructions

Assay Time ~2 hour

Properties

Form 96 well

Storage instruction Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test

reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual

for detail temperatures of the components.

Note For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Background Canine parvovirus is a contagious virus mainly affecting dogs. CPV is highly contagious and is spread

from dog to dog by direct or indirect contact with their feces. Vaccines can prevent this infection, but mortality can reach 91% in untreated cases. Treatment often involves veterinary hospitalization. Canine parvovirus often infects other mammals including foxes, wolves, cats, and skunks. Felines (cats) are also

susceptible to panleukopenia, a different strain of parvovirus.

Function Dogs that develop the disease show signs of the illness within three to ten days. The signs may include

lethargy, vomiting, fever, and diarrhea (usually bloody). Generally, the first sign of CPV is lethargy. Secondary signs are loss of weight and appetite or diarrhea followed by vomiting. Diarrhea and vomiting result in dehydration that upsets the electrolyte balance and this may affect the dog critically. Secondary infections occur as a result of the weakened immune system. Because the normal intestinal lining is also compromised, blood and protein leak into the intestines, leading to anemia and loss of

www.arigobio.com arigo.nuts about antibodies 1/2

protein, and endotoxins escape into the bloodstream, causing endotoxemia. Dogs have a distinctive odor in the later stages of the infection. The white blood cell level falls, further weakening the dog. Any or all of these factors can lead to shock and death. Younger animals have worse survival rates.