

ARG10922 anti-SV40 small t antigen antibody [pab 280]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [pab 280] recognizes SV40 small t antigen
Tested Reactivity	Virus
Tested Application	IHC-P, IP, WB
Host	Mouse
Clonality	Monoclonal
Clone	pab 280
Isotype	IgG2b
Target Name	SV40 small t antigen
Species	Virus
Immunogen	A purified Pab 419-small T complex.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	IHC-P	Assay-dependent
	IP	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified by affinity chromatography.
Buffer	PBS and 0.02% Sodium azide.
Preservative	0.02% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.