

ARG11179 anti-Nipah Virus nucleocapsid protein antibody [SQab30390]

Package: 100 µg
Store at: -20°C

Summary


Product Description	Recombinant Human Monoclonal antibody [SQab30390] recognizes Nipah Virus nucleocapsid protein
Tested Reactivity	Virus
Tested Application	ELISA
Host	Human
Clonality	Monoclonal
Clone	SQab30390
Isotype	IgG
Target Name	Nipah Virus nucleocapsid protein
Species	Virus
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein A
Purity	> 90% (SDS-PAGE)
Buffer	PBS (pH 7.4) and 0.03% Proclin 300
Preservative	0.03% Proclin 300
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG11179 anti-Nipah Virus nucleocapsid protein antibody [SQab30390] ELISA image

ELISA: The plate was coated with Nipah Virus nucleocapsid protein or BSA control. Proteins were stained with serially diluted ARG11179 anti-Nipah Virus nucleocapsid protein antibody [SQab30390].