

## ARG20620 Mouse IgG2b Isotype Control antibody [A-1] (APC-Cyanine 5.5)

Package: 50 tests  
Store at: 4°C

### Summary

Product Description	APC-Cyanine 5.5-conjugated Mouse Monoclonal antibody [A-1] as a negative control antibody for Mouse IgG2b
Tested Application	ELISA, FACS, FLISA, IHC-Fr
Specificity	Chicken IgA
Host	Mouse
Clonality	Monoclonal
Clone	A-1
Isotype	IgG2b, kappa
Target Name	Mouse IgG2b
Target Ig	Chicken IgA
Conjugation	APC-Cyanine 5.5

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	10 µl/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.