

ARG21662 Goat anti-Human IgD antibody (Texas Red), pre-adsorbed

Package: 500 µg
Store at: 4°C

Summary

Product Description	Texas Red-conjugated Goat Polyclonal antibody recognizes Human IgD
Tested Reactivity	Hu
Tested Application	FACS, FLISA, IHC-Fr, IHC-P
Specificity	The antibody reacts with the heavy chain of Human IgD. The antibody is pre-adsorbed with Human IgG, IgM and IgA, so the antibody may not react with Human IgG, IgM and IgA, but may react with IgD from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgD
Species	Human
Conjugation	Texas Red

Application Instructions

Pre Adsorbed	Human IgG, IgM and IgA.										
Application table	<table> <tr> <th>Application</th><th>Dilution</th></tr> <tr> <td>FACS</td><td>Assay-dependent</td></tr> <tr> <td>FLISA</td><td>1:100 - 1:400</td></tr> <tr> <td>IHC-Fr</td><td>Assay-dependent</td></tr> <tr> <td>IHC-P</td><td>Assay-dependent</td></tr> </table>	Application	Dilution	FACS	Assay-dependent	FLISA	1:100 - 1:400	IHC-Fr	Assay-dependent	IHC-P	Assay-dependent
Application	Dilution										
FACS	Assay-dependent										
FLISA	1:100 - 1:400										
IHC-Fr	Assay-dependent										
IHC-P	Assay-dependent										
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.										

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation

Highlight

Related news:
[Comprehensive anti-Human secondary antibodies](#)

Images

