

ARG21677

Goat anti-Human Kappa Light Chain antibody, F(ab')₂ fragment, pre-adsorbed

Package: 250 µg

Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes Human Kappa Light Chain
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, WB
Specificity	The antibody reacts with Human kappa light chains. The antibody is pre-adsorbed with Human lambda light chains, so the antibody may not react with Human lambda light chains, but may react with kappa light chains from other species.
Host	Goat
Clonality	Polyclonal
Isotype	F(ab') ₂ IgG
Target Name	Kappa Light Chain
Species	Human
Conjugation	Un-conjugated

Application Instructions

Pre Adsorbed	Human lambda light chains.	
Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties


Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw

cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.

Images



ARG21677 Goat anti-Human Kappa Light Chain antibody, F(ab')₂ fragment (pre-adsorbed) ELISA image

ELISA: The plate was coated with purified Human IgGκ, IgMκ, IgAκ, IgGλ, IgMλ, and IgAλ. Immunoglobulins were detected with serially diluted ARG21677 Goat anti-Human Kappa Light Chain antibody, F(ab')₂ fragment (pre-adsorbed) followed by [ARG23857](#) Swine anti-Goat IgG (H+L) antibody (HRP) (pre-adsorbed).