

ARG21761 Rabbit anti-Human IgG (H+L) antibody (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Rabbit Polyclonal antibody recognizes Human IgG (H+L)
Tested Reactivity	Hu
Tested Application	ELISA, FACS, FLISA, ICC/IF, WB
Specificity	Reacts with the heavy and light chains of human IgG and the light chains of human IgM and IgA
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Human
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:2000 - 1:4000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Highlight

Related news:

[Comprehensive anti-Human secondary antibodies](#)