

## ARG21763 Rabbit anti-Sheep IgG (H+L) antibody (TRITC)

Package: 500 µg  
Store at: 4°C

### Summary

Product Description	TRITC-conjugated Rabbit Polyclonal antibody recognizes Sheep IgG (H+L)
Tested Reactivity	Sheep
Tested Application	ELISA, FLISA, ICC/IF, IHC-Fr, WB
Specificity	Reacts with the heavy and light chains of sheep IgG
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Sheep
Conjugation	TRITC

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FLISA	1:100 - 1:400
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.