

## ARG21776 Pig anti-Goat IgG (H+L) antibody (PE), pre-adsorbed

Package: 250 µg  
Store at: 4°C

### Summary

Product Description	PE-conjugated Pig Polyclonal antibody recognizes Goat IgG (H+L)
Tested Reactivity	Goat
Tested Application	ELISA, FACS, FLISA, IHC-Fr, WB
Specificity	Reacts with the heavy and light chains of goat IgG
Host	Pig
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Goat
Conjugation	PE

### Application Instructions

Pre Adsorbed	Human, Rat, and Mouse serum proteins.	
Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	< 1 µg/ml
	IHC-Fr	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

**Note**

For laboratory research only, not for drug, diagnostic or other use.