

ARG21784 Mouse anti-Chicken IgM antibody [M-4]

Package: 250 µg
Store at: -20°C

Summary


Product Description	Mouse Monoclonal antibody [M-4] recognizes Chicken IgM
Tested Reactivity	Chk
Tested Application	FACS, IP
Host	Mouse
Clonality	Monoclonal
Clone	M-4
Isotype	IgM, kappa
Target Name	IgM
Species	Chicken
Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	FACS	< 1 ug/10 ⁶ cells
	IP	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21784 Mouse anti-Chicken IgM antibody [M-4] FACS image

Flow Cytometry: Chicken peripheral blood lymphocytes stained with ARG21784 Mouse anti-Chicken IgM antibody [M-4] and [ARG21153](#) anti-CD3 antibody [CT-3] (PE) followed by [ARG21518](#) Goat anti-Mouse IgM antibody (FITC) (pre-adsorbed).