

ARG21819 Mouse anti-Human IgG1 (Hinge) antibody [4E3]

Package: 250 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [4E3] recognizes Human IgG1 (Hinge)
Tested Reactivity	Hu
Tested Application	ELISA, FACS, FLISA, IHC-P, Puri, WB
Host	Mouse
Clonality	Monoclonal
Clone	4E3
Isotype	IgG1, kappa
Target Name	IgG1 (Hinge)
Species	Human
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	< 5 ug/ml
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-P	Assay-dependent
	Puri	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.