

ARG21825 Mouse anti-Human IgG1 (Fc) antibody [HP6001] (AP)

Package: 500 µl
Store at: 4°C

Summary


Product Description	AP-conjugated Mouse Monoclonal antibody [HP6001] recognizes Human IgG1 (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, FACS, FLISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	HP6001
Isotype	IgG2b, kappa
Target Name	IgG1 (Fc)
Species	Human
Immunogen	Human IgG1 myeloma protein
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:2000 - 1:8000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21825 Mouse anti-Human IgG1 (Fc) antibody [HP6001] (AP)
ELISA image

ELISA: The plate was coated with purified Human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted ARG21825 Mouse anti-Human IgG1 (Fc) antibody [HP6001] (AP).