

ARG21828 Mouse anti-Human IgG2 (Fc) antibody [31-7-4] (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated Mouse Monoclonal antibody [31-7-4] recognizes Human IgG2 (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, WB
Host	Mouse
Clonality	Monoclonal
Clone	31-7-4
Isotype	IgG1, kappa
Target Name	IgG2 (Fc)
Species	Human
Immunogen	Human IgG2 myeloma protein
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:500 - 1:1000
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	WB	Assay-dependent
	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Application Note		

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Images

