

ARG21829 Mouse anti-Human IgG2 (Fc) antibody [31-7-4] (PE)

Package: 100 µg
Store at: 4°C

Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [31-7-4] recognizes Human IgG2 (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, WB
Host	Mouse
Clonality	Monoclonal
Clone	31-7-4
Isotype	IgG1, kappa
Target Name	IgG2 (Fc)
Species	Human
Immunogen	Human IgG2 myeloma protein
Conjugation	PE

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	< 1 µg/ml
	ICC/IF	Assay-dependent
	WB	Assay-dependent
	Application Note <p>* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.</p>	

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.