

## ARG21951 Mouse anti-Human IgG2 (Fc) antibody [31-7-4]

Package: 250 µg  
Store at: -20°C

### Summary


Product Description	Mouse Monoclonal antibody [31-7-4] recognizes Human IgG2 (Fc)
Tested Reactivity	Hu
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, WB
Host	Mouse
Clonality	Monoclonal
Clone	31-7-4
Isotype	IgG1, kappa
Target Name	IgG2 (Fc)
Species	Human
Immunogen	Human IgG2 myeloma protein
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	WB	Assay-dependent
	Application Note <p>* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.</p>	

### Properties

Form	Liquid
Buffer	BBS (pH 8.2)
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21951 Mouse anti-Human IgG2 (Fc) antibody [31-7-4] ELISA image

ELISA: The plate was coated with purified Human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted ARG21951 Mouse anti-Human IgG2 (Fc) antibody [31-7-4] followed by [ARG21886](#) Goat anti-Mouse IgG1 antibody (HRP) (pre-adsorbed).