

ARG22700 Goat anti-Cat IgG (Fc) antibody (FITC)

Package: 500 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Goat Polyclonal antibody recognizes Cat IgG (Fc) This antibody recognizes feline IgG (Fc) and shows no cross-reactivity with other feline immunoglobulin classes in immunoelectrophoresis.
Tested Reactivity	Cat, Rb
Tested Application	FACS, IHC-Fr
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (Fc)
Species	Cat
Immunogen	Purified feline IgG (Fc) fragment.
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
	IHC-Fr	1:200 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.09% Sodium azide and 0.2% BSA.
Preservative	0.09% Sodium azide
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.