

ARG22970 Mouse anti-Rat IgA (Heavy Chain) antibody [MARA-1] (FITC)

Package: 250 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Mouse Monoclonal antibody [MARA-1] recognizes Rat IgA (Heavy Chain) Mouse anti Rat IgA Heavy Chain antibody, clone MARA-1 recognizes the alpha heavy chain of rat immunoglobulin. Mouse anti Rat IgA Heavy Chain antibody, clone MARA-1 shows no cross-reactivity with other rat immunoglobulin classes.
Tested Reactivity	Rat
Tested Application	FACS
Host	Mouse
Clonality	Monoclonal
Clone	MARA-1
Isotype	IgG1
Target Name	IgA (Heavy Chain)
Species	Rat
Immunogen	Purified IR1060 IgA Rat myeloma protein.
Conjugation	FITC

Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified by affinity chromatography.
Buffer	PBS, 0.1% Sodium azide and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.