

ARG23183
anti-APAAP Complex antibody [WR31]

Package: 500 µl

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [WR31] recognizes APAAP Complex APAAP complexes are designed for use as detection reagents in immunohistochemistry procedures. They are used in combination with a suitable bridging antibody and appropriate visualization system.
Tested Reactivity	Other
Tested Application	IHC-Fr, IHC-P, IHC-R
Host	Mouse
Clonality	Monoclonal
Clone	WR31
Isotype	IgG1
Target Name	APAAP Complex
Species	Others
Immunogen	Bovine alkaline phosphatase.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	IHC-Fr	1:30 - 1:50
	IHC-P	1:30 - 1:50
	IHC-R	1:30 - 1:50

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Unpurified.
Buffer	RPMI1640 medium and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.