

ARG23246 anti-MCM2 antibody [CRCT2.1]

Package: 50 µg

Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [CRCT2.1] recognizes MCM2 Mouse anti Human MCM2 antibody, clone CRCT2. 1 recognizes human DNA replication licensing factor MCM2, also known as Mcm2, Nuclear protein BM28 or mini chromosome maintenance protein-2. Mcm-2 is a 904 amino acid ~125kDa nuclear protein involved in the control of DNA replication. Mouse anti Human MCM2 antibody, clone CRCT2. 1 has been used successfully for the detection of human MCM2 in ovarian adenocarcinoma by immunohistochemistry on formalin fixed, paraffin embedded tissues (Gakiopoulou et al. 2007).
Tested Reactivity	Hu
Tested Application	IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	CRCT2.1
Isotype	IgG1
Target Name	MCM2
Species	Human
Immunogen	Bacterially expressed Human MCM2.
Conjugation	Un-conjugated
Alternate Names	CCNL1; Nuclear protein BM28; BM28; DNA replication licensing factor MCM2; cdc19; Minichromosome maintenance protein 2 homolog; MITOTIN; CDCL1; D3S3194; EC 3.6.4.12

Application Instructions

Application table	Application	Dilution
	IHC-P	1:200
	WB	Assay-dependent
Application Note	IHC-P: Antigen Retrieval: Boil tissue section in Sodium citrate buffer (pH 6.0). * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties


Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Concentration	1 mg/ml

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	MCM2
Gene Full Name	minichromosome maintenance complex component 2
Background	The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex (pre_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein forms a complex with MCM4, 6, and 7, and has been shown to regulate the helicase activity of the complex. This protein is phosphorylated, and thus regulated by, protein kinases CDC2 and CDC7. Multiple alternatively spliced transcript variants have been found, but the full-length nature of some variants has not been defined. [provided by RefSeq, Oct 2012]
Function	Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells. The active ATPase sites in the MCM2-7 ring are formed through the interaction surfaces of two neighboring subunits such that a critical structure of a conserved arginine finger motif is provided in trans relative to the ATP-binding site of the Walker A box of the adjacent subunit. The six ATPase active sites, however, are likely to contribute differentially to the complex helicase activity. Required for the entry in S phase and for cell division. [UniProt]
Calculated Mw	102 kDa
PTM	Phosphorylated on Ser-108 by ATR in proliferating cells. Ser-108 proliferation is increased by genotoxic agents. Ser-40 is mediated by the CDC7-DBF4 and CDC7-DBF4B complexes, while Ser-53 phosphorylation is only mediated by the CDC7-DBF4 complex. Phosphorylation by the CDC7-DBF4 complex during G1/S phase is required for the initiation of DNA replication. [UniProt]

Images



ARG23246 anti-MCM2 antibody [CRCT2.1] IHC-P image

Immunohistochemistry: Paraffin-embedded Human tonsil stained with ARG23246 anti-MCM2 antibody [CRCT2.1]. Antigen Retrieval: Boil tissue section in Sodium citrate buffer (pH 6.0).