

ARG23753
Mouse anti-Goat IgG (Fc) antibody [SB115h] (PE)Package: 50 µg
Store at: 4°C

Summary


Product Description	PE-conjugated Mouse Monoclonal antibody [SB115h] recognizes Goat IgG (Fc)
Tested Reactivity	Goat, Sheep
Tested Application	FLISA
Host	Mouse
Clonality	Monoclonal
Clone	SB115h
Isotype	IgG2b
Target Name	IgG (Fc)
Species	Goat
Conjugation	PE

Application Instructions

Application table	Application	Dilution
	FLISA	< 1 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG23753 Mouse anti-Goat IgG (Fc) antibody [SB115h] (PE) FLISA image

FLISA: The plate was coated with purified goat IgG Fc and IgG Fab. Immunoglobulins were detected with serially diluted ARG23753 Mouse anti-Goat IgG (Fc) antibody [SB115h] (PE).