

ARG23859

Donkey anti-Rat IgG (H+L) antibody (FITC), pre-adsorbed

Package: 500 µg

Store at: 4°C

Summary

Product Description	FITC-conjugated Donkey Polyclonal antibody recognizes Rat IgG (H+L)
Tested Reactivity	Rat
Tested Application	ELISA, FLISA, IHC-P, WB
Specificity	Reacts with the heavy and light chains of rat IgG and the light chains of rat IgM
Host	Donkey
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Rat
Conjugation	FITC

Application Instructions

Pre Adsorbed	Mouse serum proteins.	
Application table	Application	Dilution
	ELISA	Assay-dependent
	FLISA	1:200 - 1:400
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.