

ARG23861

Donkey anti-Rat IgG (H+L) antibody (HRP), pre-adsorbed

Package: 500 µl

Store at: 4°C

Summary

Product Description	HRP-conjugated Donkey Polyclonal antibody recognizes Rat IgG (H+L)
Tested Reactivity	Rat
Tested Application	ELISA, FLISA, IHC-P, WB
Specificity	Reacts with the heavy and light chains of rat IgG and the light chains of rat IgM
Host	Donkey
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Rat
Conjugation	HRP

Application Instructions

Pre Adsorbed	Mouse serum proteins.	
Application table	Application	Dilution
	ELISA	1:4000 - 1:8000
	FLISA	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4) and 50% Glycerol.
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.