

ARG23924 Mouse anti-Monkey IgG (H+L) antibody [SB108a] (PE)

Package: 50 µg
Store at: 4°C

Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [SB108a] recognizes Monkey IgG (H+L)
Tested Reactivity	Mk, R. Mk
Species Does Not React With	Hu, Ms, Rat, Bov, Cat, Chk, Dnk, Goat, Gpig, Hm, Hrs, Pig, Rb, Sheep
Tested Application	ELISA, FACS, FLISA, WB
Specificity	This antibody reacts with Rhesus and cynomolgus IgG. Minimal reactivity to human and rabbit immunoglobulins, goat IgG, and mouse, rat, hamster, guinea pig, sheep, donkey, bovine, horse, porcine, feline, and chicken serum.
Host	Mouse
Clonality	Monoclonal
Clone	SB108a
Isotype	IgG1
Target Name	IgG (H+L)
Immunogen	Rhesus monkey IgG.
Conjugation	PE

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	< 1 µg/ml
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.