

ARG24034 Goat anti-Human IgE antibody (HRP)

Package: 500 µg
Store at: 4°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Human IgE
Tested Reactivity	Hu, R. Mk
Tested Application	ELISA, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgE
Species	Human
Immunogen	Native Human IgE.
Conjugation	HRP

Application Instructions

Application table	Application	Dilution
	ELISA	1:10000 - 1:100000
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS, 0.05% Pro-Clean 400 and 0.2% BSA.
Preservative	0.05% Pro-Clean 400
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.