

ARG24090 anti-Salmonid Ig antibody

Package: 50 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Salmonid Ig This antibody recognizes salmonid immunoglobulin, binding to Ig of the Salmo, Oncorhynchus and Salvelinus genera. It may be used to detect increased levels of immunoglobulins in fish serum, indicative of parasitic or viral infection.
Tested Reactivity	Fsh
Tested Application	ELISA, ICC/IF, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Salmonid Ig
Species	Fish
Immunogen	Rainbow trout and Atlantic salmon immunoglobulins.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	1:5000 - 1:10000
	ICC/IF	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Serum
Buffer	Serum and < 0.1% Sodium azide.
Preservative	< 0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.