

ARG30284 TIM1 / KIM1 ELISA Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10621	anti-TIM1 / KIM1 antibody [KIM70]	Mouse mAb	Hu	ELISA, WB	50 µg
ARG10622	anti-TIM1 / KIM1 antibody [KIM75]	Mouse mAb	Hu	ELISA, WB	50 µg

Summary

Target Name	TIM1 / KIM1
Alternate Names	TIM1 / KIM1 ELISA antibody; TIM1 / KIM1 antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	ELISA Antibody Duo for TIM1 / KIM1
Highlight	Related Product: anti-TIM1 / KIM1 antibody;