

## ARG30286 Transferrin ELISA Antibody Duo

Package: 1 pair  
Store at: -20°C

### Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10448	anti-Transferrin antibody [2A2]	Mouse mAb	Hu	ELISA	125 µg
ARG10447	anti-Transferrin antibody [1C10]	Mouse mAb	Hu	ELISA	125 µg

### Summary

Target Name	Transferrin
Alternate Names	Transferrin ELISA antibody; Transferrin antibody

### Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Gene Full Name	ELISA Antibody Duo for Transferrin
Highlight	Related Product: <a href="#">anti-Transferrin antibody;</a>