

ARG30290 Osteocalcin ELISA Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10528	anti-Osteocalcin antibody [2H9]	Mouse mAb	Hu	ELISA	125 µg
ARG10529	anti-Osteocalcin antibody [6F9]	Mouse mAb	Hu	ELISA	125 µg

Summary

Target Name	Osteocalcin
Alternate Names	Osteocalcin ELISA antibody; Osteocalcin antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	ELISA Antibody Duo for Osteocalcin
Highlight	Related Product: anti-Osteocalcin antibody;