

ARG30291
serum Amyloid A ELISA Antibody DuoPackage: 1 pair
Store at: -20°C**Component**

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10639	anti-serum Amyloid A antibody [SAA19]	Mouse mAb	Cat	ELISA, WB	50 µg
ARG10640	anti-serum Amyloid A antibody [VSA34]	Mouse mAb	Cat	ELISA, WB	50 µg

Summary

Target Name	serum Amyloid A
Alternate Names	serum Amyloid A ELISA antibody; serum Amyloid A antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	ELISA Antibody Duo for serum Amyloid A
Highlight	Related Product: anti-serum Amyloid A antibody;