

ARG42307 Mouse anti-Human IgD antibody [IA6-2] (APC)

Package: 50 tests

Store at: 4°C

Summary

Product Description	APC-conjugated Mouse Monoclonal antibody [IA6-2] recognizes Human IgD
Tested Reactivity	Hu
Tested Application	FACS
Specificity	The mouse monoclonal antibody IA6-2 recognizes human immunoglobulin D.
Host	Mouse
Clonality	Monoclonal
Clone	IA6-2
Isotype	IgG2a, kappa
Target Name	IgD
Species	Human
Immunogen	Human IgD.
Conjugation	APC

Application Instructions

Application table	Application	Dilution
	FACS	10 µl / 100 µl of whole blood or 10 ⁶ cells
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Buffer	PBS and 15 mM Sodium azide.
Preservative	15 mM Sodium azide
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.