

## ARG42995 anti-C3orf38 antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes C3orf38
Tested Reactivity	Hu
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	C3orf38
Species	Human
Immunogen	Recombinant protein of Human C3orf38.
Conjugation	Un-conjugated
Alternate Names	Uncharacterized protein C3orf38

### Application Instructions

Application table	Application	Dilution
	WB	1:1000 - 1:5000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	K562	
Observed Size	~ 38 kDa	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	50 mM Tris-Glycine (pH 7.4), 150 mM NaCl, 0.01% Sodium azide, 40% Glycerol and 0.05% BSA.
Preservative	0.01% Sodium azide
Stabilizer	40% Glycerol and 0.05% BSA
Concentration	Batch dependent
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	C3orf38
Gene Full Name	chromosome 3 open reading frame 38
Function	May be involved in apoptosis regulation. [UniProt]
Calculated Mw	38 kDa

Images



ARG42995 anti-C3orf38 antibody WB image

Western blot: K562 cell lysate stained with ARG42995 anti-C3orf38 antibody at 1:1000 dilution.