

## ARG43486 anti-PDIA4 / ERp72 antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes PDIA4 / ERp72.
Tested Reactivity	Hu, Ms, Rat
Tested Application	ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	PDIA4 / ERp72
Species	Human
Immunogen	Synthetic peptide derived from human PDIA4 / ERp72
Conjugation	Un-conjugated
Alternate Names	ERP70; ERP72; ERp-72

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	IHC-P	1:50 - 1:200
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% sodium azide.
Preservative	0.02% sodium azide.
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation


Gene Symbol	PDIA4
Gene Full Name	protein disulfide isomerase family A, member 4

### Images




ARG43486 anti-PDIA4 / ERp72 antibody ICC/IF image

Immunofluorescence: HeLa stained with ARG43486 anti-PDIA4 / ERp72 antibody at 1:50 dilution.




ARG43486 anti-PDIA4 / ERp72 antibody WB image

Western blot: HepG2 and MCF7 stained with ARG43486 anti-PDIA4 / ERp72 antibody at 1:1000 dilution.




ARG43486 anti-PDIA4 / ERp72 antibody WB image

Western blot: Rat liver stained with ARG43486 anti-PDIA4 / ERp72 antibody at 1:1000 dilution.




ARG43486 anti-PDIA4 / ERp72 antibody IHC-P image

Immunohistochemistry: Mouse kidney stained with ARG43486 anti-PDIA4 / ERp72 antibody at 1:50 dilution.



ARG43486 anti-PDIA4 / ERp72 antibody ICC/IF image

Immunofluorescence: NIH/3T3 stained with ARG43486 anti-PDIA4 / ERp72 antibody at 1:50 dilution.



ARG43486 anti-PDIA4 / ERp72 antibody WB image

Western blot: Mouse liver stained with ARG43486 anti-PDIA4 / ERp72 antibody at 1:1000 dilution.