

ARG43994 anti-PSMC4 antibody

Package: 50 µg
Store at: -20°C

Summary


Product Description	Rabbit Polyclonal antibody recognizes PSMC4
Tested Reactivity	Hu, Ms, Rat
Tested Application	ELISA, FACS, ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	PSMC4
Species	Human
Immunogen	Human PSMC4 recombinant protein
Conjugation	Un-conjugated
Alternate Names	PSMC4; Proteasome 26S Subunit, ATPase 4; TBP-7; TBP7; Tat-Binding Protein 7; MIP224; RPT3; S6; Proteasome (Prosome, Macropain) 26S Subunit, ATPase, 4; 26S Proteasome AAA-ATPase Subunit RPT3; 26S Proteasome Regulatory Subunit 6B; MB67-Interacting Protein; Protease 26S Subunit 6; MGC13687; MGC23214; MGC8570; 26S Protease Regulatory Subunit 6B; Proteasome 26S Subunit ATPase 4; MB67 Interacting Protein

Application Instructions

Application table	Application	Dilution
	ELISA	0.1-0.5 µg/ml
	FACS	1-3 µg/1x10 ⁶
	ICC/IF	5 µg/ml
	IHC-P	2-5 µg/ml
	WB	0.1-0.25 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity purified with Immunogen.
Buffer	0.9% NaCl, 0.2% Na ₂ HPO ₄ and 4% Trehalose.
Stabilizer	4% Trehalose
Concentration	0.5 mg/ml




ARG43994 anti-PSMC4 antibody IHC-P image

Immunohistochemistry: Human adenocarcinoma stained with ARG43994 anti-PSMC4 antibody at 2 µg/ml dilution.




ARG43994 anti-PSMC4 antibody WB image

Western blot: HepG2 and U251 stained with ARG43994 anti-PSMC4 antibody at 0.5 µg/mL dilution.



ARG43994 anti-PSMC4 antibody FACS image

Flow Cytometry: U20S cells stained with ARG43994 anti-PSMC4 antibody (blue) at 1 µg/1x10⁶ cells dilution.




ARG43994 anti-PSMC4 antibody IHC-P image

Immunohistochemistry: Human glioblastoma stained with ARG43994 anti-PSMC4 antibody at 2 µg/ml dilution.

ARG43994 anti-PSMC4 antibody WB image

Western blot: Rat liver and Rat brain stained with ARG43994 anti-PSMC4 antibody at 0.5 µg/mL dilution.



ARG43994 anti-PSMC4 antibody WB image

Western blot: Mouse liver and Mouse brain stained with ARG43994 anti-PSMC4 antibody at 0.5 µg/mL dilution.

