

ARG44187 anti-MEGF9 antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes MEGF9
Tested Reactivity	Hu
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	MEGF9
Species	Human
Immunogen	Recombinant protein of Human MEGF9
Conjugation	Un-conjugated
Alternate Names	MEGF9; Multiple EGF Like Domains 9; EGFL5; Multiple Epidermal Growth Factor-Like Domains Protein 9; Epidermal Growth Factor-Like Protein 5; EGF-Like-Domain, Multiple 5; Multiple EGF-Like Domains Protein 9; Multiple EGF-Like-Domains 9; EGF-Like Protein 5; KIAA0818

Application Instructions

Application table	Application	Dilution
	WB	0.25-0.5 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.05% Sodium azide and 5% BSA.
Preservative	0.05% Sodium azide
Stabilizer	5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	MEGF9
Gene Full Name	Multiple Egf Like Domains 9
Background	Predicted to be involved in several processes, including animal organ morphogenesis; cell migration; and substrate adhesion-dependent cell spreading. Predicted to be integral component of membrane. Predicted to be active in basement membrane.
Calculated Mw	63 kDa
PTM	Disulfide bond, Glycoprotein. [Uniprot]
Cellular Localization	Membrane. [Uniprot]

Images



ARG44187 anti-MEGF9 antibody WB image

Western blot: HeLa, A549 and RT4 stained with ARG44187 anti-MEGF9 antibody at 0.5 µg/mL dilution.