

ARG45016 anti-IgG1 antibody [RM117]

Package: 50 µg
Store at: -20°C

Summary


Product Description	Rabbit monoclonal antibody [RM117] recognizes Human IgG1
Tested Reactivity	Hu
Tested Application	ELISA, FACS, ICC/IF, IHC-P
Host	Rabbit
Clonality	Monoclonal
Clone	RM117
Isotype	IgG
Target Name	IgG1
Species	Human
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	50 - 200 ng/well (for Capture)
	FACS	2 - 8 µg/ml
	ICC/IF	0.5 - 2 µg/mL
	IHC-P	0.5 - 2 µg/mL
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS with 50% Glycerol, 1% BSA and 0.09% sodium azide
Preservative	0.09% sodium azide
Stabilizer	50% Glycerol, 1% BSA and 0.09%
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG45016 anti-IgG1 antibody [RM117] IHC-P image

Immunohistochemistry: Human Tonsil stained with ARG45016 anti-IgG1 antibody [RM117].