

ARG45033 anti-25-OH Vitamin D3 antibody antibody [RM428]

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit monoclonal antibody [RM428] recognizes 25-OH Vitamin D3/RM3 or 25-OH Vitamin D3/ RMH04 Complexes
Tested Reactivity	Other
Tested Application	ELISA
Host	Rabbit
Clonality	Monoclonal
Clone	RM428
Isotype	IgG
Target Name	25-OH Vitamin D3
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	0.5 µg/mL - 2 µg/mL

Application Note