

ARG45037
anti-5-Carboxylcytosine (5-caC) antibody [RM462]

Package: 50 µg

Store at: -20°C

Summary

Product Description	Rabbit monoclonal antibody [RM462] recognizes 5-Carboxylcytosine (5-caC) Rabbit Monoclonal Antibody.
Tested Reactivity	All
Tested Application	ELISA
Host	Rabbit
Clonality	Monoclonal
Clone	RM462
Isotype	IgG
Target Name	5-Carboxylcytosine (5-caC)
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	0.2 µg/mL - 2 µg/mL
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS with 50% Glycerol, 1% BSA and 0.09% sodium azide
Preservative	0.09% sodium azide
Stabilizer	50% Glycerol, 1% BSA and 0.09%
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.