

ARG45854 anti-Zebrafish Pcdh15b antibody

Package: 200 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Pcdh15b
Tested Reactivity	Zfsh
Tested Application	IHC-P
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Pcdh15b
Species	Zebrafish
Conjugation	Un-conjugated
Alternate Names	PCDH15B; Protocadherin Related 15B;

Application Instructions

Application table	Application	Dilution
	IHC-P	2-5 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


Properties

Form	Liquid
Purification	Affinity chromatography purified
Buffer	0.02M PBS, 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% glycerol
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	PCDH15B
-------------	---------

Images



ARG45854 anti-Zebrafish Pcdh15b antibody IHC-P image

Immunohistochemistry: Zebrafish brain stained with ARG45854 anti-Zebrafish Pcdh15b antibody at 2 µg/ml dilution.