

ARG54548 anti-V5 tag antibody (Agarose)

Package: 50 µg
Store at: 4°C

Summary

Product Description	Agarose beads coupled Goat Polyclonal antibody recognizes V5 tag
Tested Reactivity	Other
Tested Application	IP
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	V5 tag
Species	Others
Conjugation	Agarose

Application Instructions

Application table	Application	Dilution
	IP	Use 15-25ul of gel slurry per 0.1-1mg of protein extract or lysate
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Buffer	PBS (pH 7.2) and 0.09% sodium azide.
Preservative	0.09% sodium azide
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.