

ARG54555 anti-VSV-G tag antibody

Package: 50 µg
Store at: 4°C

Summary

Product Description	FITC-conjugated Goat Polyclonal antibody recognizes VSV-G tag
Tested Reactivity	Other
Tested Application	ICC/IF, IHC-P
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	VSV-G tag
Species	Others
Immunogen	VSV-G (YTDIEMNRLGK) coupled to KLH
Conjugation	FITC

Application Instructions

Application Note	IHC, ICC/IF: 0.5 - 5 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Buffer	PBS (pH 7.2), 0.2% BSA and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Stabilizer	0.2% BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody
---------------	-------------------------------