

ARG55268 anti-DNA antibody [ET844.3.1]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [ET844.3.1] recognizes DNA
Tested Reactivity	All
Tested Application	ELISA
Specificity	This antibody recognize single- or double-stranded Eucaryotic DNA and E. coli DNA.
Host	Mouse
Clonality	Monoclonal
Clone	ET844.3.1
Isotype	IgM
Target Name	DNA
Species	Bovine
Immunogen	Bovine thymus DNA.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	PEG-purified
Buffer	PBS (7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Gene Regulation antibody
---------------	--------------------------