

**ARG56540**  
**anti-GST antibody (FITC)**Package: 250 µl  
Store at: 4°C**Summary**

Product Description	Rabbit Polyclonal antibody recognizes GST
Tested Reactivity	Other
Tested Application	ICC/IF, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	GST
Immunogen	Purified glutathione S-transferase (GST) from <i>S. japonicum</i> .
Conjugation	FITC

**Application Instructions**

Application table	Application	Dilution
	ICC/IF	1:200 - 1:400 (1.0 - 0.5 µg/ml)
	WB	1:200 - 1:400 (1.0 - 0.5 µg/ml)
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

**Properties**

Form	Liquid
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.