

ARG62366 anti-B Cell antibody [B381 (IPO10)]

Package: 100 µl
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [B381 (IPO10)] recognizes B Cell
Tested Reactivity	Hu
Tested Application	IHC-Fr, IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	B381 (IPO10)
Isotype	IgG
Target Name	B Cell
Species	Human
Immunogen	Human B cell line Daudi
Conjugation	Un-conjugated

Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

Properties

Form	Liquid
Purification	Protein A purified
Buffer	PBS (pH 7.2), 0.2% BSA and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Stabilizer	0.2% BSA
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Immune System antibody; Signaling Transduction antibody
---------------	---