

## ARG62372 anti-Blood Group antigen B antibody [B460]

Package: 100 µl

Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [B460] recognizes Blood Group antigen B
Tested Reactivity	Hu
Tested Application	IHC-Fr, IHC-P
Specificity	stains erythrocytes and vascular epithelium of blood group B controls.
Host	Mouse
Clonality	Monoclonal
Clone	B460
Isotype	IgM
Target Name	Blood Group antigen B
Species	Human
Immunogen	Blood Group B erythrocytes
Conjugation	Un-conjugated

### Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
Positive Control	Human Tissue (Blood Group B or AB).

### Properties

Form	Liquid
Buffer	1X PBS buffer with < 0.1% sodium azide.
Preservative	< 0.1% sodium azide.
Concentration	2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Cell Biology and Cellular Response antibody; Immune System antibody
---------------	---