

## ARG62542 anti-Medroxyprogesterone Acetate antibody [MP01(1F5)]

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rat Monoclonal antibody [MP01 (1F5)] recognizes Medroxyprogesterone Acetate
Tested Reactivity	Other
Tested Application	ELISA, RIA
Host	Rat
Clonality	Monoclonal
Clone	MP01 (1F5)
Isotype	IgG1
Target Name	Medroxyprogesterone Acetate
Species	Others
Immunogen	Medroxyprogesterone Acetate (MPA) -3-carboxymethyl oxime BSA.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	RIA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified Antibody
Buffer	1X PBS and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

