

## ARG63535 anti-SMUG1 antibody

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes SMUG1
Tested Reactivity	Hu
Tested Application	IHC-P, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	SMUG1
Species	Human
Immunogen	PQAFLGSIHEPA-C
Conjugation	Un-conjugated
Alternate Names	Single-strand selective monofunctional uracil DNA glycosylase; UNG3; HMUDG; EC 3.2.2.-; FDG

### Application Instructions

Application table	Application	Dilution
	IHC-P	Assay - dependent
	WB	0.1 - 0.3 µg/ml
Application Note	WB: Recommend incubate at RT for 1h. IHC-P: Antigen Retrieval: Steam tissue section in Citrate buffer (pH 6.0). * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	


### Properties

Form	Liquid
Purification	Purified from goat serum by antigen affinity chromatography.
Buffer	Tris saline (pH 7.3), 0.02% Sodium azide and 0.5% BSA.
Preservative	0.02% Sodium azide
Stabilizer	0.5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation


Database links	<a href="#">GeneID: 23583 Human</a> <a href="#">Swiss-port # Q53HV7 Human</a>
Background	This gene encodes a protein that participates in base excision repair by removing uracil from single- and double-stranded DNA. Many alternatively spliced transcript variants exist for this gene; the full-length nature is known for some but not all of the variants. [provided by RefSeq, Aug 2011]
Research Area	Gene Regulation antibody
Calculated Mw	30 kDa

## Images



ARG63535 anti-SMUG1 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human cortex tissue. Antigen Retrieval: Steam tissue section in Citrate buffer (pH 6.0). The tissue section was stained with ARG63535 anti-SMUG1 antibody at 5 µg/ml dilution followed by AP-staining.



ARG63535 anti-SMUG1 antibody WB image

Western blot: 35 µg of MOLT-4 cell lysate (in RIPA buffer) stained with ARG63535 anti-SMUG1 antibody at 0.1 µg/ml dilution and incubated at RT for 1 hour.