

## ARG64892 anti-FEM1A / APRAP antibody

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes FEM1A / APRAP
Tested Reactivity	Hu
Predict Reactivity	Cow, Dog
Tested Application	WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	FEM1A / APRAP
Species	Human
Immunogen	C-GFIPEDLEAFIELH
Conjugation	Un-conjugated
Alternate Names	EPRAP; Protein fem-1 homolog A; Prostaglandin E receptor 4-associated protein; FEM1-alpha; FEM1a

### Application Instructions

Application table	Application	Dilution
	WB	1 - 3 µg/ml
Application Note	WB: Recommend incubate at RT for 1h. * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties


Form	Liquid
Purification	Purified from goat serum by antigen affinity chromatography.
Buffer	Tris saline (pH 7.3), 0.02% Sodium azide and 0.5% BSA.
Preservative	0.02% Sodium azide
Stabilizer	0.5% BSA
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

---

Database links	<a href="#">GeneID: 55527 Human</a> <a href="#">Swiss-port # Q9BSK4 Human</a>
Gene Symbol	FEM1A
Gene Full Name	fem-1 homolog a (C. elegans)
Function	Probable component of an E3 ubiquitin-protein ligase complex, in which it may act as a substrate recognition subunit (By similarity). May participate in antiinflammatory signaling via its interaction with PTGER4. [UniProt]
Research Area	Cell Biology and Cellular Response antibody; Cell Death antibody
Calculated Mw	74 kDa

## Images



ARG64892 anti-FEM1A / APRAP antibody WB image

Western Blot: Human Heart lysate (35 µg protein in RIPA buffer) stained with ARG64892 anti-FEM1A / APRAP antibody at 1 µg/ml dilution.