

ARG65645 Donkey anti-Goat IgG antibody (TRITC)

Package: 1 ml
Store at: 4°C

Summary

Product Description	TRITC-conjugated Donkey Polyclonal antibody recognizes Goat IgG
Tested Reactivity	Goat
Tested Application	ELISA, ICC/IF, IHC-Fr
Specificity	Tetramethylrhodamine isothiocyanate-conjugated IgG fraction of polyclonal Donkey anti-Serum to Goat IgG, heavy and light chains.
Host	Donkey
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG
Species	Goat
Target Ig	Goat IgG
Conjugation	TRITC

Application Instructions

Cross Reactivity Note	Caprine. Inter-species cross-reactivity is a normal feature of antibodies to immunoglobulins, since Ig of different species frequently share antigenic determinants. Cross-reactivity of this antiSerum has been tested in double radial immunodiffusion with the following results: Bovine ++ Guinea Pig ± Mouse ± Rat ± Cat ± Hamster - Monkey + Sheep ++ Chicken - Horse - pigeon - Swine - Dog ± Human + Rabbit ± Turkey -								
Application table	<table> <tr> <th>Application</th><th>Dilution</th></tr> <tr> <td>ELISA</td><td>Assay-dependent</td></tr> <tr> <td>ICC/IF</td><td>Assay-dependent</td></tr> <tr> <td>IHC-Fr</td><td>Assay-dependent</td></tr> </table>	Application	Dilution	ELISA	Assay-dependent	ICC/IF	Assay-dependent	IHC-Fr	Assay-dependent
Application	Dilution								
ELISA	Assay-dependent								
ICC/IF	Assay-dependent								
IHC-Fr	Assay-dependent								
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.								

Properties

Form	Liquid
Buffer	PBS (pH 7,2)
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Research Area	Immune System antibody
---------------	------------------------